

HydROC Fiber and Hydrolysis Unit

One instrument for both Fiber and Hydrolysis Easy to use Fiber bag system All samples separated



LP 1023 www.liquidline.se

specifications subject to change without notice



The sample is weighed into an OPSIS LiquidLINE Fiberbag.



Sealing of bag with Thermosealer ensures that sample remains inside the bag during the whole process



Boiling according to Weende Fiber method.



Efficient water rinsing, using the HydROC drain valves.



Removal of Filterbag for ashing and weighing.

OPSIS AB, Box 244 SE-244 02 Furulund Sweden Telephone +46 46 72 25 00 Telefax +46 46 72 25 01 E-mail info@opsis.se www.liquidline.se





LP1023 2021 02

